


Center-open Vise Grip GTW-50

A sample, regardless of the size, is always gripped in the center of the vise as the grip open and shut at the center.





Lever type

Features

- The grips always hold a test sample at the center of the vise.
- The jaw faces have a vertical V-groove that can firmly hold a round shaft.
- Easy to operate dial knob type and strong grip lever type are available.

Note

- Some samples may not correctly be gripped depending on material and shape.
- Information within this document is subject to change without notice.

84464 Waldkraiburg Germany

Phone: int. + 49 / (0)8638 / 9410-0 int. + 49 / (0)8638 / 4825 Fax:

e-mail: info@hans-schmidt.com Internet: http://www.hans-schmidt.com



Specifications

	Туре	Capacity	Jaw face width	Weight
GTW-50R	Dial knob	500 N (50 kgf)	50 mm	3 kg
GTW-50L	Lever	1000 N (100 kgf)	50 mm	3 kg

Dimensions

